

CSCA 5642 Final Project

Convolutional Neural Networks for Photoelectric-
Compton Decomposition in Dual-Energy CT Tissue
Classification

Karl Schmidt

Table of Contents

1. Introduction
2. Methodology
3. Evaluation & Analysis
4. Conclusion

INTRODUCTION

Introduction

- Goal: Create a DL-spectral CT model to accurately “recover” three-tissue images.
- Models:
 - 11 convolutional layers (Unet256) and
 - 15 convolutional layers (UNet512)

METHODOLOGY

Methodology


Dataset

Dual-Energy CT Scan (Spectral CT)

- Low Energy 50kVp
- High Energy 80kVp




Dataset – for 50 & 80 kVp

- 1000x512x512 True Tissue Values
- 1000x256x1024 kVp Transmission





Methodology

Exploratory Data Analysis



Methodology


Exploratory Data Analysis



Methodology

Data Preparation


Convert the transmission data into Attenuation values.



$$A = \log_{10} \left(\frac{I_0}{I} \right) = \epsilon c l$$


Use Filtered BackProjection (FBP) to reconstruct kVp Images from Attenuation values $[\mu_{\text{low}}, \mu_{\text{high}}]$.

Filtered Back Projection



Methodology

Models



UNet256

- 2 encoding and decoding levels
- Encodes up to 256 pixels
- Faster to train but likely less expressive

UNet512

- 3 encoding and decoding levels
- Encodes up to 512 pixels
- Potential to recorder greater image complexity.

```
class ConvBlock(nn.Module):
    def __init__(
        self,
        in_channels,
        out_channels,
        kernel_size=3,
        padding=1
    ):
        super().__init__()
        self.block = nn.Sequential(
            nn.Conv2d(in_channels, out_channels, kernel_size=kernel_size, padding=padding),
            nn.ReLU(inplace=True),
            nn.Conv2d(out_channels, out_channels, kernel_size=kernel_size, padding=padding),
            nn.ReLU(inplace=True),
        )

    def forward(self, x):
        return self.block(x)
```

EVALUATION & ANALYSIS

Evaluation & Analysis

Hyperparameter Optimization

Tuned Hyperparameters

- **Filter Size:** 3, 4, 5
- **Stride:** 2, 3
- **Padding:** 1, 2
- **Learning Rate:** 1e-4, 1e-3, 1e-2
- **Batch Size:** 4, 8
- **Epochs:** 3, 5

Best Configurations

- UNet256: [Val Loss: 0.0519]
 - Filter Size: 4
 - Stride: 3
 - Padding: 2
 - Learning Rate: 0.001
 - Batch Size: 4
- UNet512: [Val Loss: 0.0380]
 - Filter Size: 5
 - Stride: 3
 - Padding: 1
 - Learning Rate: 0.001
 - Batch Size: 4

Evaluation & Analysis

Final Training Results

Epoch	Train Loss	Val Loss	MAE (Adipose, Fibro, Calc)
1	0.2254	0.0456	(0.00420, 0.00260, 0.00023)
2	0.0478	0.0366	(0.00469, 0.00469, 0.00010)
3	0.0365	0.0331	(0.00161, 0.00175, 0.00005)
4	0.0324	0.0310	(0.00216, 0.00260, 0.00004)
5	0.0304	0.0292	(0.00260, 0.00278, 0.00003)

Loss Curves

- Train Loss: Decreased from 0.2254 to 0.0304
- Val Loss: Dropped from 0.0456 to 0.0292 over 5 epochs.

Epoch	Predicted (%)			Ground Truth (%)		
	Adipose	Fibro	Calc	Adipose	Fibro	Calc
1	39.35	22.46	0.06	39.67	22.35	0.04
2	40.12	21.89	0.03	39.67	22.35	0.04
3	39.51	22.52	0.05	39.67	22.35	0.04
4	39.83	22.14	0.05	39.67	22.35	0.04
5	39.42	22.61	0.04	39.67	22.35	0.04

Evaluation & Analysis

Interpretation

Key Takeaways

- UNet512 Outperforms UNet256
- Model converges quickly
- High tissue composition accuracy

Challenges & Considerations

- Class imbalance
- Low error doesn't guarantee clinical interpretability

Conclusion

UNet512 demonstrates effective learning of dual-energy spectral features for multi-tissue segmentation in noiseless 2D breast CT.

CONCLUSIONS

Conclusion

Results Summary

- UNet 512 achieved the best performance
 - Final validation loss: 0.0292
 - Tissue MAE: as low as <0.003
 - Tissue compositions closely matched ground truth
- Demonstrated robust performance despite class imbalance and small calcification signals

Conclusion

Key Contributions

- Integrated Beer-Lambert attenuation law into training design
- Used physically meaningful $[\mu_{\text{low}}, \mu_{\text{high}}]$ features for segmentation
- Trained and evaluated fully end-to-end U-Net models on simulated AAPM data

Conclusion

Future Directions

- Add data augmentation or class balancing for rare tissue types (e.g., calcification)
- Extend model to noisy or real clinical data
- Explore multi-task learning (e.g., decomposition + classification)